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Emergency Stretcher

使用说明书
Users' Manual

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1 Product Feather

Emergency stretcher adopts materials like aluminum alloy, with the feather of light structure, smart, durable, easy to be adjusted from high to low position and easy to load and unload from the ambulance.

2 Main application and scope

Used for medical institutions to transport and move patients.


3 Structher

The stretcher is consisted of stretcher surface, frame, guardrail, wheel, safety belt and optional accessories. Optional accessories include: mattresses, IV pole. Passive products.

4 Product model division description and basic parameters (Refer to table 1, table 2)

4.1 Product model: WSX-AC-1, WSX-AC-2, WSX-AC-3, WSX-AC-5

4.2 Product model division description

4.3 Basic parameters

Table 1 Basic parameters

Unit:mm

Specific ation and model	Structure	Bed surface size		Guardrail size		Upper wheel height (Stretcher deployment status)
		Length	Width	Length	Height from bed surface	
WSX-AC-1	Stretcher surface and frame are detachable; The height of the frame can be adjusted in three gears; Stretcher surface with adjustable leg rest for lifting; The stretcher has two folding operating handles at the rear end	1950±50	550 ⁺⁵⁰ ₋₂₀	≥500	175±25	≤750
WSX-AC-2	The stretcher body can be converted into a chair posture; There is no folding operating handle at the back of the stretcher					
WSX-AC-3	There is no folding operating handle at the back of the stretcher					
WSX-AC-5	There is one folding operating handle at the back of the stretcher					

Table 2 Basic parameters

Unit:mm

Project		WSX-AC-1	WSX-AC-2	WSX-AC-3	WSX-AC-5
Detachable stretcher connection height	The height from the lowest contact surface of the detachable stretcher and its frame to the lying surface of the stretcher	≤ 300	-----	-----	-----
Regulation	Backrest length	≥ 600			
	Lift angle of backrest part	$\geq 45^\circ$	$\geq 60^\circ$	$\geq 75^\circ$	$\geq 75^\circ$
	Lifting angle of leg rest	$\geq 10^\circ$	-----	-----	-----
Frame size	After folding the stretcher, the length of the frame should not exceed the length of the stretcher surface	100			
	After folding the stretcher, the width of the frame should not exceed the width of the stretcher surface	50			

5 Loading capacity:The Max capacity 180kgs.

6 Installation, adjustment, and use

6.1 After unpacking, push the stretcher into the ambulance carriage and select a suitable position to fix the front and rear fixing devices with bolts on the floor of the car carriage. If the ambulance is equipped with a stretcher platform, the front and rear fixing devices can be installed on the stretcher

platform according to this method.

⚠ Warning: When pushing the wheels of the stretcher onto the floor of the ambulance carriage, it is important to ensure that there is a ground clearance (3–5 cm) between the front wheels of the stretcher, otherwise the front legs of the stretcher will not be able to open, leading to injury to the patient!

⚠ Attention: During installation, please increase the distance between the front and rear fixing by 5mm appropriately, otherwise there may be a possibility that the stretcher cannot be fixed in place after a period of use!

6.2 Patient gets on the stretcher

6.2.1 Support one side of the guardrail and lower the other side of the guardrail(Figure 1).

6.2.2 Lift the patient onto the stretcher, ensuring that their head is on the backrest side of the stretcher(Figure 2).

6.2.3 Raise the guardrail and fasten the seat belt(Figure 3).

6.2.4 If you need to adjust the angle of the back panel, one must hold the back panel with one hand and press the adjustment handle with the other hand. Release the handle to automatically lock in the gear position.

Figure 1

Figure 2

Figure 3

6.3 Stretcher gets on the ambulance

6.3.1 The main operator pushes the stretcher, and make the upper wheels of the stretcher fully reach the ambulance floor or above the stretcher platform(Figure 5),the deputy operator cooperates with the main operator to push and lift the stretcher at the back board to prevent the stretcher from tipping and the patient from falling(Figure 6).

6.3.2 Getting on ambulance process

- WSX-AC-1 model:After the front leg of the stretcher touches the ambulance bumper, pull the left operating handle of the stretcher to fold the front leg of the stretcher, and simultaneously push the stretcher until the rear leg contacts the ambulance bumper. Pull the right operating handle to fold the rear leg of the stretcher, and continue to push the stretcher until it is fully pushed into the ambulance(Figure 7).
- WSX-AC-2,WSX-AC-5 model:After the front leg of the stretcher touches the bumper of the ambulance, pull the stretcher operating handle and simultaneously push the stretcher until it is fully pushed into the ambulance(Figure 7).
- WSX-AC-3 model:After the front leg trigger rod of the

stretcher touches the bumper of the ambulance, push the stretcher until it is fully pushed into the ambulance (Figure 7).

⚠ Attention: When loading the stretcher onto the ambulance, the main operator must push/lift the stretcher to prevent it from slipping or tipping.

Figure 5

Figure 6

Figure 7

6.3.3 After the stretcher is pushed into the ambulance carriage, adjust the position of the stretcher slightly to make the front and rear fixing devices are in place.

⚠ Attention: After pushing the stretcher into the ambulance, it is necessary to check and ensure that the front and rear fixing devices of the stretcher are fixed in place, otherwise the stretcher may detach from the ambulance during operation.

6.4 Stretcher gets off the ambulance

Getting off Ambulance process

- WSX-AC-1, WSX-AC-3, WSX-AC-5 model: The main operator unlocks the rear fixing device, slowly pulling the stretcher out and lifting it slightly upwards until the rear legs of the

stretcher are fully open and touching the ground; Continue to slowly pull out the stretcher, at which point the deputy operator must support the front of the stretcher until the front legs of the stretcher are fully open and touching the ground. Continue to pull out the stretcher and the upper wheels are off the ambulance floor.

- WSX-AC-2 model: The main operator unlock the rear fixing device and slowly pull the stretcher out while lifting it slightly upwards until the front and rear legs of the stretcher are fully open and touching the ground. Continue to pull out the stretcher and the upper wheels are off the ambulance floor.

⚠ Attention: Do not touch the stretcher operating handle during the stretcher is getting off the ambulance.

⚠ Warning: When the stretcher getting off the ambulance, before the upper wheels on the stretcher detach from the ambulance floor, it is important to ensure that the front and rear legs of the stretcher are fully open, otherwise there is a risk of injury to the patient and damage to the stretcher!

6.5 Split and install of stretcher surface and frame

6.5.1 Split of stretcher surface and frame

- Pull the unlocking handle and pull the stretcher surface

backwards (Figure 8)

- The front and rear two people lift the stretcher surface and detach it from the frame.


Figure 8

6.5.2 Installation of stretcher surface and frame

The front and rear two people lift the stretcher surface onto the stretcher frame and push in to lock it.

⚠ Attention: After installation, make sure that the stretcher surface and frame are securely locked, otherwise the stretcher surface and frame may detach!

6.6 Frame height adjustment (WSX-AC-1 model is applicable)

Pull the operating handle of the stretcher leg rest to adjust the height of the frame (Figure 9).

Pull the operating handle at the backrest end of the stretcher to adjust the height of the frame (Figure 10).


Figure 9

Figure 10

⚠ Attention:When adjusting the height of the frame, both front and rear person must support the stretcher to avoid the danger of the stretcher surface falling!

⚠ Attention:It is not recommended to adjust the frame when there is patient on the stretcher.

6.7 Leg rest adjustment (WSX-AC-1 model is applicable)

Lift/pull up and secure the leg rest

6.8 The stretcher is converted from a bed position to a chair position (WSX-AC-2 model is applicable)

Hold the lock release mechanism with one hand and press it down, while holding the upper wheel end with the other hand to lift it up. Adjust the angle of the back plate as needed (Figure 11).


Figure 11

⚠ Attention:Regardless of whether the stretcher is in a bed

or chair position, it must be ensured that the unlocking mechanism is locked in place after adjustment, otherwise it may damage the patient or the stretcher!

⚠ Attention: Only when there are no patient on the stretcher can displacement adjustment be performed!

⚠ Warning: During displacement adjustment, all parts of the human body should not be close to or under the stretcher, otherwise there is a danger of crushing!

7. Maintenance and disinfection

7.1 Use facial cloth or soft brush, clear the stretcher with little water and neutral detergent and then dry it with a clean face cloth.

7.2 Use neutral disinfectant spray or wiping disinfection, never use Iodine or phenolic disinfectants. After disinfection by spray, the stretcher shall be dried immediately.

7.3 Inspect weekly , tighten the screws of various parts.

7.4 Fill lubricating oil to rotating parts every month.

8. Notice

8.1 To ensure the convenient and safe transportation of patient with stretcher, it is recommend that the using unit allocate dedicated stretcher operators who should receive training.

Untrained person are not allowed to operate this stretcher!

8.2 At least two people must operate this stretcher before and after. When rescuing the injured, it is necessary to pull up the guardrail, fasten the safety belt, and strictly prohibit rough or irregular operations, otherwise there is a possibility of falling the patient and damaging the stretcher! (Figure 12)

8.3 The position, quantity, and operation method of the operating handles of different models of stretchers vary. At the beginning of use, it is necessary to conduct operational drills using an unmanned stretcher. Personnel who have not receive the above learning and drills are not allowed to operate this stretcher!

8.4 On the basis of the stretcher being able to smoothly get on and off the ambulance, it is necessary to ensure that the height of the ambulance floor is slightly higher than the height of the boarding wheel from the ground (3-5 cm), otherwise the front legs of the stretcher may not be fully opened when getting off the ambulance, and there is a possibility of injury to the patient and damage to the stretcher!

8.5 When the stretcher is fixed on the ambulance, it is necessary to check whether the front and rear fixation are intact and effective, and ensure that the stretcher is fixed

in place, otherwise there is a possibility of the stretcher detaching from the ambulance while it is moving.

8.6 When pulling down the stretcher from the ambulance, the deputy operator must support the front of the stretcher until the front legs of the stretcher are fully open and the upper wheels on the stretcher are off the ambulance floor, otherwise there is a possibility of falling and damaging the stretcher!\

8.7 It is strictly prohibited to push on bumpy and uneven ground, touch the folding mechanism of the front and rear support legs of the stretcher, and violently collide with ambulances or other objects (Figure 13).

8.8 When getting on or off the ambulance, please choose a flat ground and strictly prohibit getting on or off the ambulance on slopes or bumpy roads!

8.9 This product is only suitable for transporting patients. It is strictly prohibit to overload or ship other items (Figure 14)

8.10 Do not rinse (Figure 15)

8.11 Regular inspections should be conducted to ensure that bolts in all parts are not loose and lubricating oil should be added to all moving parts.

8.12 If any parts are damaged, please contact the dealer or our company's after-sales department. Using non company parts will pose a safety hazard when using stretchers!

8.13 Unauthorized modification of this stretcher will pose a safety hazard when using it!

8.14 The damage to stretcher components, operation methods, and usage environment can all pose safety hazards for stretcher use. Once the front and rear legs of the stretcher cannot be fully opened, the operators in front and behind the stretcher must firmly support the stretcher to prevent it from falling and causing injury to the patient.

Figure 12

Figure 13


Figure 14





Figure 15

9. Contraindication: None

10. Graphics and symbols

Graphics and symbols	Meaning	Graphics and symbols	Meaning
	Batch number		Place upwards during transportation or storage
	Identification of product		Afraid of rain

 Warning	When used incorrectly, there is a possibility of personal injury or death		Limit of stacking layers
 Attention	When used incorrectly, there is a possibility of equipment damage and personal injury or death		Fragile items, handle with care during handling
	Manufacturer		Authorised representative in the European Community

11. Transportation and Storage

When shipping, general transportation or according to the order contract can be used. During transportation, severe collision, heavy pressure, and exposure to sunlight and rain should be avoided; The packaged emergency stretcher should be stored indoors in a well ventilated and non corrosive atmosphere. Avoid mixing with toxic and harmful substances.

12. Contact information and after-sale service

The warranty period of this product is one year, and the following situations are not covered by the warranty:

- ① Damage or failure caused by misused.
- ② Damage or failure caused by failure to perform maintenance according to this manual.
- ③ Damage or failure caused by unforeseeable natural disasters.

④ The wear and tear of vulnerable parts, such as deformation of mattresses, wear and tear of wheels and fenders.

During the warranty period, our company provides free repair parts, and after the warranty period, we are also responsible for product repair/parts replacement and charge corresponding fees. Not responsible for other quality assurance responsibilities beyond warranty.

Accessories List

Name	Quantity	Replacement cycle and method
Safety belt	2	Purchase a new safety belt
Mattress	1	Purchase a new mattress
Fixing device	1	Loosen the four screws on the fixing device and remove them. Replace them with a new fixing device and tighten the screws again.
Certificate of conformity	1	/
Instructions	1	/


Manufacturer

Name: Jiangyin Everise Medical Equipment Co.,Ltd

Add: 1001 Chengyang Road,Jiangyin City,Jiangsu,China

Tel: +86 400 158 1358

www.Everisemedical.com

Authorised representative in the European Community

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Share Info Consultant Service LLC Repräsentanzbüro